


Biosensis® LipoFluor-P1™ TR-600-P1™ Ready-to-Dilute™, Polar Lipid Tracing Reagent

Catalog Number: Lipofluor-P1™ TR-600-P1



For research use only, not for use in clinical and diagnostic procedures.

1. Intended Use

Amphiphilic polar lipids are found in membranes and perform many biological functions. Most polar lipids in cell membranes are glycerophospholipids (GPLs) with a glycerol backbone and fatty acids. GPLs include phosphatidylcholine, ethanolamine, and serine. The polar lipids derived from sphingosin are called sphingophospholipids. The most representative neural SPL is sphingomyelin, which contains ceramide and phosphorylcholine. Dietary polar lipids alter membrane compositions and signaling and enzymatic functions because they interact with cellular membranes. Polar lipids have anti-inflammatory, anti-cancer, cholesterol-lowering, and brain-development effects.

Lipofluor-P1™TR-600-P1 localizes with polar lipids, allowing the labeling of lipid droplets and other high lipid content compartments in various cell and tissue types. Lipofluor-P1™ TR-600-P1 is compatible with live and fixed samples and can be imaged by single-photon or multi-photon microscopy. Lipofluor-P1™ TR-600-P1 is also compatible with vibrational spectroscopy as this complex is infrared and Raman active.

Tested Applications: Live cell cultures and unfixed tissues; fixed cells and tissues.

Validated systems:

Cell penetration and localization of TR-600-P1 have been confirmed in a range of cell lines, including adipocytes (3T3-L1), prostate cells (PNT2, PNT1a, LNCaP, 22RV1, and DU145), cardiomyocytes (H9c2) and neuronal cells (PC-12), and tissues, including adipose tissue (sheep and Drosophila), muscle tissue (sheep cardiac and skeletal), and brain (murine).

2. Materials Provided

Four 0.5 mL microfuge vials containing 5.5 uL/vial of 10 mM stock solution of Lipofluor-P1™ TR-600-P1 in DMSO. 34.1 µg per vial. Each vial can be diluted 100X to 2000X depending upon the final working concentration of the dye in the experiment. Typical working concentrations for Lipofluor-P1™ TR-600-P1 are 10µM-50µM.

3. Specifications

CAS Number: 1369583-75-4; FW: 620.6

MF: C₂₃H₁₂N₇O₃Re

Ex/Em: 405 nm/ 550 nm: Form: Liquid:

Working concentrations: 10µM-20µM

4. Precautions for Use

It is not recommended that detergents such as Tween20 or supplements with high lipid content, such as fetal calf serum, be used to prepare these reagents as the lipid in the solution will bind the dye, reducing the cell signal. Low solubility in aqueous solutions may cause dye precipitate if used at concentrations higher than recommended. Before performing the staining procedure for fixed or live cell imaging, please read the entire procedure and consider the safety data sheet. For laboratory use only. Not fully tested. Not for drug, household, human, or veterinary uses.

5. Reagent Preparation

Note: Lipofluor-P1™ TR-600-P1 undiluted stock solution vials should be stored at 2-8°C protected from light for up to 12 months from receipt. Diluted working solutions are short-lived. Dilute to working solution just before use. Excess solutions do not store well. See “Storage conditions” below for more information.

To create a 10 µM, 5.445 mL working solution, follow these steps for each 0.5 mL vial:


- Take a 0.5 mL vial containing 5.5 µL of LipoFluor™ TR-600-P1 in DMSO liquid.
- Dilute the DMSO liquid with 49.5 µL of an aqueous solution (such as PBS or cell media) to create a 10x dilution.
- Transfer this mixture to a larger vial and dilute it further with 5.445 mL of the aqueous solution (100x dilution) to make a total of 5.5 mL of Lipofluor-P1™ TR-600-P1 solution at 10 µM. This operating solution is enough to generate solutions for more than 24 wells per vial with 200 µL per well.
- If you need to make a specific solution concentration, you can use a basic chemistry formula to calculate how much of the original stock solution to use. The formula is $C_1V_1=C_2V_2$, where C_1 is the concentration of the stock solution, V_1 is the volume of the stock solution needed, C_2 is the desired concentration of the final solution, and V_2 is the final volume of the solution. Thus, one will solve for V_1 's volume. For example, if you need to make 1 milliliter of a 50 µM solution using a 10 mM stock solution, you can use an online calculator (see figure) or solve the C formula manually. In this case, C_1 is 10 mM, C_2 is 50 µM, and V_2 is 1 mL. Solving for V_1 , you get $V_1=C_2V_2/C_1$, which means adding 5.0 µL of the stock solution to 995 µL of solvent to make ~1 mL of a 50 µM solution.

6. Storage Conditions

Lipofluor-P1™ TR-600-P1 stock 10 mM DMSO solution vials should be stored, tightly sealed, and well protected from light at 2-8°C for up to 12 months from receipt. Opened stock tubes are stable if stored tightly sealed and protected from light for six months from purchase. Working or diluted solutions are unstable and should be made just before use. We have no data on frozen solutions; thus, we cannot recommend freezing the DMSO stocks or diluted solutions.

Calculating Molar Dilutions:

For example, this online source:
<https://bit.ly/3leWFbR>



4. Dilute a stock solution

Stock concentration: 10 millimolar

Desired concentration: 50 micromo

Desired volume: 1 milliliter

Required volume = 5 microliter

7. Staining Protocol For Live Cells

For Adherent Cells

1. Grow cells in 6 well-plate on coverslips with appropriate culture medium and under appropriate growth conditions
2. Grow cells to the desired confluence (70 – 80%)
3. Remove culture medium and add pre-warmed PBS, pH 7.2-7.6, or cell culture media containing 10 – 20 µM of Lipofluor-P1™ TR-600-P1 (1:500-1000 dilution of 10 mM stock solution)
4. Incubate cells for 30 minutes under appropriate growth conditions
5. Wash coverslips twice for one minute in PBS, pH 7.2-7.6
6. Mount coverslips in aqueous mounting media for imaging

Note: Glycerol-based mounting media may reduce the fluorescence intensity of Lipofluor-P1™ TR-600-P1.

For Suspended Living Cells

1. Centrifuge cell suspension to obtain cell pellet and remove the supernatant
2. Resuspend cells in pre-warmed PBS, pH 7.2-7.6 (37°C) or serum-free medium containing 10 – 20 µM of Lipofluor-P1™ TR-600-P1 (1:500-1:1000 dilution of 10 mM stock solution)
3. Incubate cells for 30 minutes under appropriate growth conditions
4. Re-pellet the cells by centrifugation and resuspend in PBS, pH 7.2-7.6, or cell culture medium
5. Cells can be prepared as a wet mounted or adhere to poly-L-lysine coated coverslips

and mounted in an aqueous mounting media for immediate imaging

For Co-Staining Live cell Experiment

1. Before co-staining, ensure the spectral profiles of the counter-staining agent and Lipofluor-P1™ TR-600-P1 can be appropriately resolved. Generally, dyes that do not excite with 405 nm excitation can be imaged alongside Lipofluor-P1™ TR-600-P1. Blue dyes such as DAPI are also compatible as they emit at a lower wavelength than Lipofluor-P1™ TR-600-P1
2. Stain cells as described above with a reduced washing step to 30 seconds following incubation
3. Stain cells with a counter-staining agent according to the manufacturer's instructions
4. Following washes, mount in an aqueous mounting media for imaging

1. Staining Protocol For Fixed Cells Cell Fixation

We recommend to fix samples in 4% paraformaldehyde for 20 minutes at room temperature. Wash samples 3 x 10 minutes in PBS, pH 7.2-7.6.

Staining Fixed Cells

1. Incubate fixed cells in PBS, pH 7.2-7.6, containing 10-20 µM TR-600-P1 for 30 minutes at room temperature.
2. For best results, provide gentle agitation by a platform rocker (or similar) at low rpm.
3. Wash coverslips twice for one minute in PBS, pH 7.2-7.6, with agitation.

4. Mount coverslips in an aqueous mounting media and image immediately for best results.

8. Staining Protocol For Tissue Sections

Unlike the conventional mitochondrial stains, paraformaldehyde-fixed, fresh, and fresh-frozen tissue sections have been successfully stained with Lipofluor-P1™ TR-600-P1. Other fixation methods have not been attempted to date.

If endogenous fluorescence is an issue in your tissue sample, quenching can assist in imaging. We recommend incubating samples in 100 mM glycine in PBS (pH to 7.4 with 1 M tris base, if required) for quenching endogenous fluorescence for 20 minutes at room temperature. Other treatments, such as UV irradiation, may also help quench endogenous fluorescence; however, avoid harsh treatments that may leach lipids from samples or interfere with lipid binding.

Sample Preparation

Tissues can be stained immediately upon collection or stored for later staining. We recommend 4% paraformaldehyde fixation or flash freezing for tissue storage. Sample preparation will depend on the tissue type and imaging platform. In general, Lipofluor-P1™ TR-600-P1 can stain tissue sections of up to 5 mm thick. Live samples can be sectioned using a sharp scalpel or knife. Fixed and frozen can also be prepared in OCT sectioned by microtome to your desired thickness.

Staining Sections

Incubate samples with 10 - 20 μM Lipofluor-P1™ TR-600-P1 in PBS, pH 7.2-7.6, or

appropriate media (1:500-1:1000 dilution of 10 mM stock solution) for 30 minutes at room temperature with gentle agitation provided by a platform rocker (or similar) at low rpm. Wash samples three times for five minutes in PBS, pH 7.2-7.6, at room temperature with agitation. Mount tissue in aqueous mounting media and image immediately for best results. Optimization of both concentration and incubation times will need to be determined for best results.

9. Fluorescent Imaging Settings


Epi-Fluorescence Microscopy

Lipofluor-P1™ TR-600-P1 can be excited by UV (~365 nm) or blue light (405 nm) sources with emissions collected using a wideband pass filter or narrowband pass filter within an emission range of 550-650 nm.

Confocal or Two-Photon Microscopy

Lipofluor-P1™ TR-600-P1 can be excited by a 400 nm steady-state laser or at 800-830 nm using a two-photon pulse laser. Ideally, image the specimen with a spectral detector set for the emission of Lipofluor-P1™ TR-600-P1, 490-670 nm ($E_{max} = 570$ nm). Alternatively, detect by using an emission filter suited for detecting FITC-based fluorophores.

Note: Time-gated imaging can be performed with this product and is ideal for samples with high levels of endogenous fluorescence. Probe emission lifetime is ~30 microseconds.



Vibrational Spectroscopy

TR-600-P1 can be detected by infrared and Raman spectroscopy methods. Sample preparation will need to be adapted appropriately for these techniques.


Sample preparation

Adherent cell culture can be grown on silicon nitride substrate, which is compatible with infrared, or calcium fluoride slides, which are compatible with Raman spectroscopy. Fixation with cold methanol is recommended for best preservation of lipids. Following staining, dehydration may be required as these techniques can be hindered by water content. (Bader CA, *et al.* Mol Biosyst 2016; 12:2064-8).

Infrared Spectroscopy

The FTIR spectrum of solid TR-600-P1 is dominated by the strong Re carbonyl stretching bands centered at $\sim 2027\text{ cm}^{-1}$ and the doublet $1915/1893\text{ cm}^{-1}$, characteristic of facial tricarbonyl complexes. This is a spectral region

where vibrational modes due to biochemical components are conveniently negligible.



Raman spectroscopy

TR-600-P1 excitation can be achieved at 785 nm with a spectra as shown below.

