

# **biosensis<sup>®</sup> Bovine A1 Beta Casein**

## **ELISA Kit Protocol**

**Catalog No: BEK-2364-1P**

*For quantitative detection of bovine A1 beta casein level in bovine milk samples only  
when used as directed.*

### **TABLE OF CONTENTS**

<b>I</b>	Experimental Overview.....	2
<b>II</b>	Materials provided.....	2
<b>III</b>	Equipment required but not supplied.....	3
<b>IV</b>	Technical hints.....	3
<b>V</b>	Storage of kit components.....	4
<b>VI</b>	Reagent preparations.....	4
<b>VII</b>	Sample Preparation and Dilution.....	5
<b>VIII</b>	Assay procedure.....	6
<b>IX</b>	Calculation of results.....	7
<b>X</b>	Sample Data.....	8

## I. Experimental Overview:

The Biosensis Bovine A1 Beta Casein ELISA Kit has been designed for the quantitative detection of bovine A1 beta casein level in bovine milk samples only when used as directed. The principal of the ELISA is a sandwich assay. Affinity-purified rabbit anti-beta casein pAbs are pre-coated on ELISA plate. A1 beta casein standards (6.25 ng/mL-400 ng/mL) and tested samples will be captured and the A1 specific chicken pAbs will be used to detect the A1 beta casein. Rabbit anti-chicken IgY HRP conjugate will be used for the signal amplification. After color development, OD450 will be measured and A1 standard curve will be established. Concentration of A1 beta casein in tested sample will be estimated based on the standard curve. In general, house milk sample need to be diluted to 1:20,000 to 1:100,000 since A1 level is high. For testing A1 contamination in A2 milk, sample can be diluted to 1:10-1:1,000.

For each plate, it is recommended that two columns will be used for standard curve establishment with the detecting range of 6.25 ng/mL to 400 ng/mL A1 and that samples be run in duplicate at least.

Assay Range: 6.25 ng/mL - 400 ng/mL; Specificity: Bovine casein A1 protein, no cross-reactivity detectable with A2 allele.

---

## II. Materials Provided

Kit Components	Unit Size/Volume	Qty in Kit
ELISA Plate coated with rabbit anti-beta Casein pAbs	96 well	1
A1 beta casein standard (40 µg/mL / 40 ng/µL)	20 µL	1
Chicken anti-A1 beta casein specific antibody (0.5 mg/mL) (2000X)	10 µL	1
Rabbit Anti-Chicken IgY HRP conjugate (1000X)	20 µL	1
BSA for Ab dilution buffer in 125 mL bottle	1 g	1
10X TBS in 125 mL bottle	100 mL	1
Tween-20 (Polyoxyethylene-Sorbitan Monolaurate), 50% solution, 1 mL tube	1 mL	1
TMBS Substrate in 30 mL bottle	6 mL	1
ELISA Record Sheet and insert instructions	--	1

**WARNING: DO NOT MIX REAGENTS SUCH AS ANTIBODY BUFFERS ETC BETWEEN DIFFERENT ASSAYS IF USING MULTIPLE KITS. USE ONLY THE BUFFERS, STANDARDS & SOLUTIONS INCLUDED FOR YOUR SPECIFIC LOT. MIXING SOLUTIONS FROM DIFFERENT TARGET KITS CAN LEAD TO POOR KIT PERFORMANCE**

All components with -20°C or 4°C are shipped on blue ice packs and will arrive cold or room temperature. See page 4 storage instructions for correct storage of kit and its components upon receipt.

---

*For Research Use Only. Not for use in Clinical & Diagnostic Procedures*

2

---

### III. Equipment Required but Not Supplied

1. Automated plate washer is desirable
2. Multi-channel or repeating pipette.
3. Clean tubes and Eppendorf tubes.
4. Plate shaker (300-500 rpm)
5. Mechanical Vortex.
6. Microplate reader with 450 nm filter.
7. Pipettors and tips
8. HCl, 1M for stopping reaction in final step

---

### IV. Technical Hints

1. Spin all kit components to get liquid to the bottom of tubes before use.
2. Duplicate well assays are recommended for both standard and sample testing. All assay steps are performed at room temperature.
3. In order to avoid marginal effects of plate incubation due to temperature differences, it is suggested that the TMB solution be brought to room temperature.
4. Recommended method for manual plate washing.
  - Emptying of wells on the plate:  
Place the plate on the palm of the hand in a position that enables easy flicking movement using the wrist. Holding the plate over a sink, quickly invert the plate, whilst accelerating the arm downward toward the sink. Abruptly stop the downward acceleration to force the liquid from the wells into the sink. When done correctly the technique should prevent liquid from getting on to the fingers or on the outside of the strip wells or plate holder. *Note: Retain the upside down position of the plate to avoid any back flow into the wells. DO NOT LET THE PLATE DRY OUT.*
  - Blotting the plate:  
Immediately blot the inverted plate by lightly tapping the plate 3-4 times on blotting paper.
  - Washing: Forcefully pipette Wash Buffer into each well with a multi-channel pipette. Empty the wells of wash buffer using technique described above. Repeat washing and flicking procedures thrice. When washing the plate, each time let the washing buffer stay in the wells for 1-2 min before rinsing with more wash buffer as directed above. *Note: **Avoid** touching the inside surface of the wells with the pipette tips.*

***Do not let the wells dry out at any time or kit performance will be poor***

---

## V. Storage of Kit Components:

The unopened kit can be stored at 4°C for up to one 1 week if the kit will be used within that week. **If the kit will not be used with 1 week of receipt, then the -20°C components should be stored in -10°C to -20°C freezer.** Expiration of the unopened kit is 6 months if stored as directed below. Once opened, use all components within 1 week. Place components in appropriate temperature locations as directed below: *Avoid repeated freeze thaw cycles on frozen items.*

ELISA Plate coated with rabbit Anti-beta Casein pAbs, 4°C  
A1 beta casein standard, 40 µg/mL (40 ng/µL), -20°C  
Chicken Anti-A1 beta casein specific pAbs, 0.5 mg/mL, -20°C  
Rabbit Anti-Chicken IgY HRP conjugate, use at 1:1,000 dilution, -20°C  
BSA for Ab dilution buffer in 125 mL bottle (need to add 100 mL TBST) 4°C  
10X TBS, 100 mL, 4°C  
Tween-20 (Polyoxyethylene-Sorbitan Monolaurate), 50%, 1 mL, 4°C  
Substrate TMBS, 4°C  
ELISA Record Sheet RT

---

## VI. Reagent Preparations

### Prepare TBS/TBST Solutions:

1. Prepare TBS: add 100 mL 10X TBS into 900 mL D.I. water, mix well; store at RT; filter through a 0.45 filter unit for lowest backgrounds.
2. Prepare TBST: Add 1 mL of the Tween-20 solution into the 1X TBS from step 1, and mix well, TBST will be used for preparing the antibody dilution buffer and for plate washing.

### Prepare Antibody/Sample Dilution Buffer:

Add 100 mL TBST solution into the 125 mL BSA bottle and mix well. Do not over mix.

### Preparation of Standards:

(a) Prepare 400 ng/mL A1 standard (tube #1): Add 10 µL of the 40 µg/mL (40 ng/µL) A1 stock protein standard into 990 µL Ab dilution buffer to make 400ng/ml A1 standard. Mix gently. This is tube #1, the high standard.

(b) 200 ng - 6.25 ng/mL standards: Label 6 additional Eppendorf tubes with (tube #2) 200 ng/mL, (tube #3) 100 ng/mL, (tube #4) 50 ng/mL, (tube #5) 25 ng/mL, (tube #6) 12.5 ng/mL, and (tube #7) 6.25 ng/mL respectively. Aliquot 0.5 mL of antibody/sample buffer into tubes 2-7. Tube 1 already has the 400 ng/mL A1 standard in it from step (a). Add

*For Research Use Only. Not for use in Clinical & Diagnostic Procedures*

4

0.5 mL from tube #1 tube into tube #2 and mix. Serially transfer 0.5 mL from 2<sup>st</sup> tube into 3<sup>nd</sup> tube and mix. Continue performing serial dilutions through tube seven. In the end, tube#1 has 400ng/ml, tube#2: 200ng/ml, #3: 100ng/ml; #4: 50ng/ml; #5: 25ng/ml; #6: 12.5ng/ml; #7: 6.25ng/ml. Store as directed below. Use sterile technique.

*Note: Standard solution should be prepared no more than 2 hours prior to the experiment. The working standard solution may be stored at 4°C for up to 3 hours. Avoid repeated freeze-thaw cycles.*

#### **Preparation of Detection Antibody:**

Prepare (10 mL per plate) A1 specific detecting pAb (Chicken IgY) in Ab dilution buffer at ~ 1:2,000 dilution, ie, add 5 µL A1 Detecting Abs into 10 mL Ab dilution buffer just before use. Do not store. See procedural instructions.

#### **Preparation of Anti-Chicken HRP Secondary:**

Prepare Rabbit Anti-Chicken IgY HRP conjugate (10 mL per plate) solution: dilute the HRP conjugate to ~ 1:1,000 in Ab dilution buffer (10 µL HRP conjugate into 10 mL Ab dilution buffer, mix well) just prior to use. Do not store. See procedural instructions.

---

## **VII. Sample Preparation and Dilution**

#### **Preparation of Samples:**

Prepare the tested milk samples by diluting milk in Ab dilution buffer accordingly. In general, A1 milk (household organic milk) can be diluted to 1:20,000-1:100,000 fold. For testing possible A1 contaminated A2 milk, you may dilute that milk 1:10-1:1,000 fold. Samples can be stored at 4°C for up to 3 hours; do not freeze.

---

## VIII. Assay Procedure

We recommend that standard solution dilutions and each sample are plated in duplicate. Read entire protocol before beginning. Prepare all buffer solutions ahead of time as directed except for the standards, actual samples, and secondary antibody dilutions. These should be made just before use.

1. Remove ELISA plate from storage and bring to room temperature. Remove the seal from ELISA plate, empty the solution and wash 3X with TBST, 1 minute each. Then add 100  $\mu$ L/well of standard or above prepared sample in duplicate and make sure to record each well location on the ELISA Record Sheet
2. Mix the plate on an ELISA Plate shaker at 400 rpm for 2 hour at RT.
3. Wash 3X with TBST, 3 minutes each.
4. While washing the plate, prepare (10 mL per plate) A1 specific detecting pAb (Chicken IgY) in Ab dilution buffer at  $\sim$  1:2,000 dilution, i.e., add 5  $\mu$ L A1 Detecting Abs into 10 mL Ab dilution buffer.
5. Add the prepared A1 detecting pAb at 100  $\mu$ L/well to the plate and shake at 400 rpm for 1 hour at RT.
6. Wash 3X with TBST, 3 minutes each.
7. While washing the plate, prepare Rabbit Anti-Chicken IgY HRP conjugate (10 mL per plate) solution: dilute the HRP conjugate to  $\sim$  1:1,000 in Ab dilution buffer (10  $\mu$ L HRP conjugate into 10 mL Ab dilution buffer, mix well). Add 100  $\mu$ L/well to the plate.
8. Mix the plate on the shaker (400 rpm) for another 1 hour at RT.
9. Wash 5X with TBST, 3-5 minutes each.
10. Add 50  $\mu$ L/well substrate (TMBS). Let the color to develop 4-8 minutes. Note the order and arrangement of the substrate addition. One must add the stop solution in the same order for best results. Color development will be quick.
11. Stop the reaction with 50  $\mu$ L 1M HCl to each well, in the same pattern and order that was used to load the TMBS substrate. It is recommended to immediately read the plate at 450 nm with a microplate reader. Note: Color will fade over time; hence, we recommend plate to be read within 30 min after adding the stop solution.

*CAUTION: Bubbles in the wells will cause inaccurate readings. Ensure that all bubbles are removed prior to taking the absorbance reading.*

12. Record the development time on the ELISA recording sheet along with OD450 readout. Make a standard curve using the average OD450 from the two reference standards for each point. Determine the A1 concentration for the tested samples based on their OD450 and standard curve. Make sure that the final concentration is obtained by multiplying with dilution factor.

---

*For Research Use Only. Not for use in Clinical & Diagnostic Procedures*

6

---

## IX. Calculation of Results

### (a) Manual Plate Reading:

(The relative O.D.<sub>450</sub>) = (the O.D.<sub>450</sub> of each well) – (the O.D.<sub>450</sub> of Zero well).

The standard curve can be plotted as the relative O.D.<sub>450</sub> of each standard solution (Y) vs. the respective concentration of the standard solution (X). Known concentrations of the target protein are plotted on the X-axis and the corresponding O.D.<sub>450</sub> on the Y-axis. The standard curve should result in a graph that shows a direct relationship between target protein concentrations and the corresponding O.D.<sub>450</sub>. The greater the concentration of target protein in the sample, the higher the O.D.<sub>450</sub>.

*Note: if the samples measured were diluted, multiply the dilution factor to the concentrations from interpolation to obtain the concentration before dilution.*

### Determine concentration of target protein in unknown sample:

The target protein concentration of the samples can be interpolated from the standard curve. Draw a horizontal line to intersect with the standard curve. A vertical line dropped from this point intersects the X-axis at the concentration of antibody in the unknown sample.

### (b) PC Interface Plate Reading

Use appropriate software to reduce the data and generate a four parameter logistic (4-PL) curve-fit; avoid using linear regression analysis. Perform a 4-PL regression analysis to calculate the concentration of target analyte in the samples. Multiply the result by the sample dilution factor.

---

## X: Typical Standard Curve (for reference only)

Room temperature, 6 minutes


A1 Std (ng/mL)	OD1	OD2	Average OD
400	1.8402	1.8738	1.8570
200	1.3734	1.3975	1.3855
100	0.7821	0.8051	0.7936
50	0.4697	0.4745	0.4721
25	0.3110	0.3154	0.3132
12.5	0.2231	0.2315	0.2273
6.25	0.1875	0.1862	0.1869
0	0.1418	0.1406	0.1412

---

*For Research Use Only. Not for use in Clinical & Diagnostic Procedures*

7

A2 400 ng/mL	0.1418	0.1409	0.1414
A1 Milk 1:20,000	1.6533	1.6178	1.6356
A2 Milk 1:1,000	0.1057	0.1036	0.1047



Graph demonstrates the dynamic range and high specificity of the beta Casein A1 ELISA. Note that samples containing A2 allele of beta casein fail to generate a signal even at 400 ng/ml concentrations (blue data point).

Biosensis.com is your antibody and ELISA kit specialist  
Helping you make sense of your research!

We offer antibodies for research in Apoptosis, Autophagy, Cancer, Development, Diabetes & Obesity, DNA Damage, GPCRs, Growth Factors, Immunology, Matrix Molecules, Neuroscience and Stem Cells.

We have over 150 highly sensitive ELISA kits with a strip format offering the highest level of sensitivity and you only need to use what you need.

**Biosensis reagents are available online directly from [www.biosensis.com](http://www.biosensis.com)**

*For Research Use Only. Not for use in Clinical & Diagnostic Procedures*

8

Biosensis Pty Ltd · 51 West Thebarton Road, Thebarton, South Australia, 5031  
Telephone +61 08 8352 7711 · Email [sales@biosensis.com](mailto:sales@biosensis.com) · [www.biosensis.com](http://www.biosensis.com)

BEK-2364-1P/3-11-2019/LMF