



# Multi-Neurotrophin *Rapid*<sup>TM</sup> Screening ELISA Kits

For the detection of *multiple* Neurotrophins  
in your samples at one Go!

Multi-Neurotrophin Standard Curves



4 in 1!



Target	Code	Species		NGF	BDNF	NT3	NT4/5
Multi-Neurotrophin	BEK-2227	Hu	Sensitivity (pg/mL)	3	1	4	30
			Range (pg/mL)	3.9 - 250	7.8 - 500	15.6 - 1000	15.6 - 1000
Multi-Neurotrophin	BEK-2231	Ms	Sensitivity (pg/mL)	3	1	4	30
			Range (pg/mL)	3.9 - 250	7.8 - 500	15.6 - 1000	15.6 - 1000
Multi-Neurotrophin	BEK-2232	Rat	Sensitivity (pg/mL)	3	1	4	30
			Range (pg/mL)	3.9 - 250	7.8 - 500	15.6 - 1000	15.6 - 1000

## Specific detection of NGF, BDNF, NT3 and NT4/5!

- Low sensitivity! Recently published! (Lindsay SL, et al. 2016;S.C.R., 6:729)
- Detection of multiple Neurotrophins with results in **less than 4hrs!**
- Excellent for a quick screen on precious samples



www.biosensis.com