



Chicken polyclonal antibody to human NT3

Catalogue No.:	C-1513-500
Description:	FUNCTION: Seems to promote the survival of visceral and proprioceptive sensory neurons. SUBCELLULAR LOCATION: Secreted protein. TISSUE SPECIFICITY: Brain and peripheral tissues. SIMILARITY: Belongs to the NGF-beta family.
Batch No.:	See product label
Unit size:	500 uL
Antigen:	Mixture of two human NT3 peptides (90-104 and 199-214 aa).
Sequence:	TELLRQRRYNSPRV, RPKVNGCRGIDDKHWN
Other Names:	Neurotrophin-3; Neurotrophic factor; HDNF; Nerve growth factor 2; NGF-2; Ntf3; Ntf-3;
Accession:	P20783 NTF3_HUMAN;
Produced in:	Chicken
Applications:	ELISA at suggested dilution of 1:500-1:2,000. Biosensis recommends that the optimal working dilution should be determined by the end user.
Cross-reactivity:	Human, Mouse
Form:	Liquid. PBS with 0.02% Sodium Azide
Storage:	Short term storage at 2-8C for one week. At -20C as an undiluted liquid for up to 12 months.
Expiry Date:	12 months after purchase

FOR RESEARCH USE ONLY