



Chicken polyclonal antibody to Obestatin (76-98): Affinity purified

Catalogue No.:	C-1530-200
Description:	Obestatin is generated from the proteolytic cleavage of Ghrelin.
Batch No.:	See product label
Unit size:	200 uL
Antigen:	Obestatin peptide (76-98 aa) from mouse ghrelin precursor.
Sequence:	FNAPFDVGIKLSGAQYQQHGRAL
Other Names:	Appetite-regulating hormone; Growth hormone secretagogue; Growth hormone-releasing peptide; Motilin-related peptide; M46 protein; Ghrelin-27; Ghrelin-28;
Accession:	Q9EQX0 GHRL_MOUSE;
Produced in:	Chicken
Applications:	ELISA . Suggested dilution at 1:500 to 1:2,000 in coating buffer as coating antibody in competitive ELISA testing. Biosensis recommends that the optimal working dilution should be determined by the end user.
Cross-reactivity:	Mouse, Rat
Form:	Liquid. PBS with 0.02% Sodium Azide
Storage:	Short term storage at 2-8C for one week. At -20C as an undiluted liquid for up to 12 months.
Expiry Date:	12 months after purchase

FOR RESEARCH USE ONLY