

Rabbit polyclonal antibody to Histone H3 (121-136): Affinity purified

Catalogue No.:	R-1549-100
Description:	Histone H3 is a core component of nucleosome.
Batch No.:	See product label
Unit size:	100 uL
Antigen:	A peptide from the C-terminus of human Histone H3 (121-136 aa).
Sequence:	MPKDIQLARRIGERA
Other Names:	Histone H3.1; HIST1H3;
Accession:	P68431 H31_HUMAN;
Produced in:	Rabbit
Applications:	WB. Suggested dilution of 1:1,000-1:5,000. Biosensis recommends that the optimal working dilution should be determined by the end user.
Specificity:	Antibody binds to various Histone H3 proteins.
Cross-reactivity:	Human, Mouse, Rat, Bovine
Form:	Liquid. PBS with 0.02% Sodium Azide
Storage:	Short term storage at 2-8C for one week. At -20C as an undiluted liquid for up to 12 months.
Expiry Date:	12 months after purchase



Western Blot analysis of Rabbit polyclonal antibody to Histone H3 (R-1549-100) in calf thymus protein extract.

FOR RESEARCH USE ONLY