


## Rabbit antibody to NT3 (139-149): affinity purified

<b>Catalogue No.:</b>	R-175-20
<b>Description:</b>	NT3 is a member of the neurotrophin family, that controls survival and differentiation of visceral and proprioceptive sensory neurons. NT3 is closely related to both NGF and BDNF. It may be involved in the maintenance of the adult nervous system, and may affect development of neurons in the embryo when it is expressed in human placenta. NT3-deficient mice generated by gene targeting display severe movement defects of the limbs. The mature peptide of this protein is identical in all mammals examined including human, pig, rat and mouse. SUBCELLULAR LOCATION: Secreted protein. TISSUE SPECIFICITY: Brain and peripheral tissues. SIMILARITY: Belongs to the NGF-beta family.
<b>Batch No.:</b>	See product label
<b>Unit size:</b>	20 ug
<b>Antigen:</b>	A synthetic peptide (YAEHKSHRGEY) as part of human (aa: 139-149), mouse and rat NT3 protein conjugated to BSA has been used as the immunogen.
<b>Other Names:</b>	Neurotrophin-3; Neurotrophic factor; HDNF; Nerve growth factor 2; NGF-2; Ntf3; Ntf-3
<b>Accession:</b>	NT3_HUMAN  NT3_MOUSE  NT3_RAT
<b>Produced in:</b>	Rabbit
<b>Purity:</b>	Affinity purified
<b>Applications:</b>	IHC, ELISA, WB, inhibition of biological activity. A concentration of 1-10 ug/mL is recommended for IHC, ELISA, WB and inhibition of biological activity in vitro; 2-10 ug/mL (ED50) for in vivo use. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
<b>Specificity:</b>	A cross reactivity of less than 0.1% to mouse NGF, recombinant human BDNF and NT4/5 has been shown by dot blot.
<b>Cross-reactivity:</b>	This antiserum is known to react with rat, chicken and human NT3.
<b>Form:</b>	Lyophilised
<b>Reconstitution:</b>	Reconstitute in 20 uL of sterile water. Centrifuge to remove any insoluble material.
<b>Storage:</b>	After reconstitution keep aliquots at -20C for a higher stability, and at 2-8C with an appropriate antibacterial agent. Glycerol (1:1) may be added for an additional stability. Avoid repetitive freeze/thaw cycles.
<b>Expiry Date:</b>	12 months after purchase
<b>References:</b>	1. Zouh et al. (1996) Mol. Neurobio. 13(3): 185-97 2. Rush et al. (1997) Clin. Exp. Pharmacol. 24: 549-55 3. Tafreshi et al. (1998) Neurosci. 83(2): 373-80

---

FOR RESEARCH USE ONLY

## Rabbit antibody to NT3 (139-149): affinity purified



Immunohistochemical staining of Neurotrophin 3 (NT3) in rat purkinje cells of the cerebellum using rabbit antibody to mouse NT3, catalogue number R-175-20.

biosensis

---

FOR RESEARCH USE ONLY