


## Mouse monoclonal antibody to human Protocadherin beta-15 [2E10]: IgG

<b>Catalogue No.:</b>	M-902-100
<b>Description:</b>	Protocadherin beta-15 is a member of the protocadherin beta gene cluster which contains 16 genes and 3 pseudogenes each encoding 6 extracellular cadherin domains and a cytoplasmic tail that deviates from others in the cadherin superfamily. Protocadherin beta-15 is a potential calcium-dependent cell adhesion protein. It may be involved in the establishment and maintenance of specific neuronal connections in the brain. It is a single-pass type I membrane protein which contains 6 cadherin domains.
<b>Unit size:</b>	100 ug
<b>Antigen:</b>	Partial recombinant human Protocadherin beta-15 [207-306] with a GST tag.
<b>Clone:</b>	[2E10]
<b>Other Names:</b>	PCDH-beta-15; PCDHB15
<b>Accession:</b>	PCDBF_HUMAN
<b>Produced in:</b>	Mouse
<b>Purity:</b>	Protein G purified immunoglobulin
<b>Applications:</b>	This antibody is recommended for WB and ELISA. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
<b>Specificity:</b>	Specificity has been confirmed by WB and ELISA against the antigen.
<b>Cross-reactivity:</b>	Human. Other species have not been tested.
<b>Form:</b>	Lyophilised from PBS pH 7.2
<b>Reconstitution:</b>	Reconstitute in 100 uL of sterile water. Centrifuge to remove any insoluble material.
<b>Storage:</b>	After reconstitution keep aliquots at -20C for higher stability or at 2-8C with an appropriate antibacterial agent. Glycerol (1:1) may be added for additional stability. Avoid repetitive freeze/thaw cycles.
<b>Expiry Date:</b>	12 months after purchase



Western blot detection of GST tagged recombinant human Protocadherin beta-15. Note the GST tag alone is 26 kDa.

FOR RESEARCH USE ONLY