

## Rabbit polyclonal antibody to human Estrogen Receptor beta (27-50): Affinity purified

<b>Catalogue No.:</b>	R-1084-100
<b>Description:</b>	Estrogen Receptor beta (ESR2) belongs to the nuclear hormone receptor family. ESR2 binds estrogens with high affinity, similar to Estrogen Receptor alpha (ESR1). At least 8 isoforms of ESR2 are produced by alternate splicing.
<b>Batch No.:</b>	See product label
<b>Unit size:</b>	100 ug
<b>Antigen:</b>	A synthetic peptide (HGSYIPSSYVDSHHEYPAMTFYS) corresponding to a region (27-50) from human Estrogen Receptor beta. To enhance the immunological response, this peptide was coupled to carrier protein BSA.
<b>Other Names:</b>	ER-beta; Nuclear receptor subfamily 3 group A member 2; ESR2; ESTRB; NR3A2;
<b>Accession:</b>	Q92731 ESR2_HUMAN;
<b>Produced in:</b>	Rabbit
<b>Purity:</b>	Affinity purified on antigen column
<b>Applications:</b>	Western Blotting (WB). A concentration of 0.5-1.0 ug/mL is recommended for WB. Human Estrogen Receptor beta (isoform Beta-1) has a predicted length of 530 residues and MW of 59 kDa. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
<b>Specificity:</b>	The specificity of this antibody has been confirmed by WB against the antigen.
<b>Cross-reactivity:</b>	Human. The antibody is predicted to react with rat due to sequence homology.
<b>Form:</b>	Lyophilised with 5mg BSA, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg Thimerosal, 0.05mg NaN <sub>3</sub>
<b>Reconstitution:</b>	Reconstitute in 100 uL of sterile distilled water to achieve an antibody concentration of 1 mg/mL. Centrifuge to remove any insoluble material.
<b>Storage:</b>	At least 12 months after purchase at 2-8C (lyophilized formulations). After reconstitution, aliquot and store at -20C for a higher stability. Avoid freeze-thaw cycles
<b>Expiry Date:</b>	12 months after purchase

---

FOR RESEARCH USE ONLY