



Mouse monoclonal antibody to G1/S-specific cyclin-D1 [CY-D1]: IgG

Catalogue No.:	M-1196-100
Description:	G1/S-specific cyclin-D1 (CCND1) is a member of the Cyclin family which act as regulators of CDK kinases. CCND1 is essential for the control of the cell cycle at the G1/S (start) transition.
Batch No.:	See product label
Unit size:	100 ug
Antigen:	Recombinant human G1/S-specific cyclin-D1 protein
Clone:	CY-D1
Other Names:	PRAD1; BCL-1 oncogene; CCND1;
Accession:	P20248 CCND1_HUMAN;
Produced in:	Mouse
Purity:	IgG
Applications:	Immunohistochemistry (IHC) and Western Blotting (WB). A concentration of 4.0 ug/mL is recommended for WB. Human G1/S-specific cyclin-D1 has a predicted length of 295 residues and MW of 34 kDa. A concentration of 8.0 ug/mL is recommended to detect the protein in formalin fixed and paraffin embedded tissues as well as formalin/acetone fixed frozen tissues. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
Specificity:	The specificity of this antibody has been confirmed by WB and IHC against the antigen.
Cross-reactivity:	Human; mouse; rat;
Form:	Lyophilized from 1.2% sodium acetate, 2mg BSA, 0.01mg NaN3
Reconstitution:	Reconstitute in 1 mL of PBS (pH 7.4) to achieve an antibody concentration of 100 ug/mL. Centrifuge to remove any insoluble material.
Storage:	At least 12 months after purchase at 2-8C (lyophilized formulations). After reconstitution, aliquot and store at -20C for a higher stability. Avoid freeze-thaw cycles
Expiry Date:	12 months after purchase

FOR RESEARCH USE ONLY