



Human Survivin ELISA Kit (2 plates)

Catalogue No.: BEK-2121-2P

Description: The human Survivin Kit is a sandwich ELISA. The capture antibody is a polyclonal human Survivin antibody pre-coated onto the 96-well strip plates provided in the kit. Human test samples and standards of known Survivin concentration are added to these wells and allowed to complex with the bound Survivin antibody. A biotinylated human Survivin polyclonal antibody is then added. This detection antibody binds to the antigen thus completing the sandwich. After washing, an enzyme Avidin-Biotin-Peroxidase complex (ABC) is added which binds to the second antibody. The peroxidase substrate TMB is added to induce a coloured reaction product. The intensity of this coloured product is directly proportional to the concentration of Survivin present in the samples. The purpose of this kit is the in-vitro quantitative determination of human Survivin in samples such as sera, plasma, tissue lysates and cell culture supernates. This kit has been configured for research use only and is not to be used in diagnostic or clinical procedures.

Batch No.: See product labels

Other Names: Baculoviral IAP repeat-containing protein 5; Apoptosis inhibitor 4; Apoptosis inhibitor survivin; BIRC5; API4; IAP4;

Accession: O15392 BIRC5_HUMAN;

Specificity: Human Survivin

Storage: Store at 4°C

Kit components: The ELISA kit box contains 2 x 96 pre-coated strip plates, protein standards, detection reagents, substrate buffer and precise instructions.


Range: 62.5 pg/ml - 4,000 pg/ml

Sensitivity: < 2 pg/ml

Kit protocol: Please refer to our online product listing for current protocol/MSDS versions.

Typical Human Survivin ELISA Kit Standard Curve
(TMB reaction incubated at 37°C for 22 min)

Concentration (pg/ml)	0.0	62.5	125	250	500	1000	2000	4000
O.D.	0.027	0.148	0.273	0.517	0.939	1.486	1.963	2.165



This standard curve is for demonstration purposes only. A standard curve should be generated for each assay.

FOR RESEARCH USE ONLY