

Parkin control peptide (305-323)

Catalogue No.:	PE-1253-100
Batch No.:	See product label
Unit size:	100 ug
Sequence:	RILGEE QYNRYQQYGA EEC
Other Names:	Ubiquitin E3 ligase PRKN; Parkinson juvenile disease protein 2; Parkinson disease protein 2; PARK2; PRKN
Accession:	PRKN2_HUMAN
Purity:	85% or greater
Applications:	Control peptide. This peptide may be used to block binding activity of antibody to Parkin (# R-113-100). Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
Form:	Lyophilised
Reconstitution:	This peptide is soluble in PBS and water.
Storage:	After reconstitution keep at -20C. Avoid repetitive freeze/thaw cycles.

biosensis

FOR RESEARCH USE ONLY
