




# Mature BDNF/proBDNF Combo Rapid™ ELISA Kits



- Combo Rapid™ kits offer the best of both kits and save you \$\$\$, too!
- Two complete kits in one Box! All the reagents and QC controls for both analytes are included.
- Two excellent kit options for measuring BDNF isoforms independently in human blood or rodent samples.
- The inclusion of the heterophilic antibody blocker assures accurate proBDNF quantification in human blood.
- Mature BDNF ELISA assay has low cross-reactivity (~5.3%) with proBDNF, while the proBDNF ELISA has little to no cross-reactivity with mature BDNF.
- Our Combo Rapid™ kits provide the same high accuracy, speed, and reproducibility that you expect from our individual kits but at a more accessible and affordable price point for today's researchers.



## Mature BDNF Rapid ELISA Kit Catalog No: BEK-2211

For the preferential detection of Mature BDNF independently of the proBDNF form.



## proBDNF Rapid ELISA Kit Catalog No: BEK-2237

Detection of BDNF in full length. The mature form is not detected.




Mature BDNF Rapid™ ELISA Kit		proBDNF Rapid™ ELISA Kit	
<b>Target</b>	Mature BDNF	<b>Target</b>	proBDNF
<b>Code</b>	BEK-2211-1P	<b>Code</b>	BEK-2237-1P
<b>Species Reactivity</b>	Hu, Ms, Rat	<b>Species Reactivity</b>	Hu, Ms, Rat
<b>Sensitivity</b>	1-2 pg/mL	<b>Sensitivity</b>	6 pg/mL
<b>Range</b>	7.8 pg/mL - 500 pg/mL	<b>Range</b>	15.6 pg/mL - 1000 pg/mL

## Cat No. 2241-SET

Optimized for human blood samples




# Mature BDNF/proBDNF Combo Rapid™ ELISA Kits



- Combo Rapid™ kits offer the best of both kits and save you \$\$\$, too!
- Two complete kits in one Box! All the reagents and QC controls for both analytes are included.
- Two excellent kit options for measuring BDNF isoforms independently in human blood or rodent samples.
- The inclusion of the heterophilic antibody blocker assures accurate proBDNF quantification in human blood.
- Mature BDNF ELISA assay has low cross-reactivity (~5.3%) with proBDNF, while the proBDNF ELISA has little to no cross-reactivity with mature BDNF.
- Our Combo Rapid™ kits provide the same high accuracy, speed, and reproducibility that you expect from our individual kits but at a more accessible and affordable price point for today's researchers.



## Mature BDNF Rapid ELISA Kit Catalog No: BEK-2211

For the preferential detection of Mature BDNF independently of the proBDNF form.



## proBDNF Rapid ELISA Kit Catalog No: BEK-2217

Detection of BDNF in full length. The mature form is not detected.



Mature BDNF Rapid™ ELISA Kit		proBDNF Rapid™ ELISA Kit	
<b>Target</b>	Mature BDNF	<b>Target</b>	proBDNF
<b>Code</b>	BEK-2211-1P	<b>Code</b>	BEK-2217-1P
<b>Species Reactivity</b>	Hu, Ms, Rat	<b>Species Reactivity</b>	Hu, Ms, Rat
<b>Sensitivity</b>	1-2 pg/mL	<b>Sensitivity</b>	6 pg/mL
<b>Range</b>	7.8 pg/mL - 500 pg/mL	<b>Range</b>	15.6 pg/mL - 1000 pg/mL

## Cat No. 2240-SET

Optimized for lysates and rodent samples

