



NGF *Rapid*TM ELISA Kits (Human, Mouse and Rat)



NGF *Rapid*TM ELISA Kits


- Biosensis has redeveloped **faster** and more **sensitive** NGF ELISA Kits for Human, Mouse and Rat NGF.
- Now you can assay with high sensitivity and accuracy and get your results in **less than 4hrs!**
- All of our ELISA Kits are available in **1-Plate** or **2-Plate** (96/192 well) formats.
- No cross reactivity with BDNF, NT3, NT4/5, GDNF or VEGF.

Standard Curves




- Representative standard curves of our NGF *Rapid*TM ELISA kits

Linearity of Dilutions



- NGF concentrations in human plasma samples using our Human NGF *Rapid*TM ELISA Kit (BEK-2212-2P).

Spike and Recovery in Rat Brain



- Spike and recovery experiment with NGF using our Rat NGF *Rapid*TM ELISA Kit (BEK-2214-2P).

Biosensis NGF *Rapid*TM ELISA Kits

Target	Code (2 Plate)	Code (1 Plate)	Species reactivity	Sensitivity	Range
NGF <i>Rapid</i> TM	BEK-2212-2P	BEK-2212-1P	Hu	2 pg/mL	3.9 pg/mL - 250 pg/mL
NGF <i>Rapid</i> TM	BEK-2213-2P	BEK-2213-1P	Ms	8 pg/mL	3.9 pg/mL - 250 pg/mL
NGF <i>Rapid</i> TM	BEK-2214-2P	BEK-2214-1P	Rat	8 pg/mL	3.9 pg/mL - 250 pg/mL