


Human Granulocyte colony-stimulating factor, G-CSF/Pluripoietin ELISA Kit (2 plates)

Catalogue No.:	BEK-2019-2P
Description:	The human Granulocyte Colony-Stimulating factor (G-CSF) Kit is a sandwich ELISA. The capture antibody is a polyclonal human G-CSF antibody pre-coated onto the 96-well strip plates provided in the kit. Human test samples and standards of known G-CSF concentration are added to these wells and allowed to complex with the bound G-CSF antibody. A biotinylated human G-CSF polyclonal antibody is then added. This detection antibody binds to the antigen thus completing the sandwich. After washing, an enzyme Avidin-Biotin-Peroxidase complex (ABC) is added which binds to the second antibody. The peroxidase substrate TMB is added to induce a coloured reaction product. The intensity of this coloured product is directly proportional to the concentration of G-CSF present in the samples. The purpose of this kit is the in-vitro quantitative determination of human G-CSF in samples such as sera, plasma, tissue lysates and cell culture supernates. This kit has been configured for research use only and is not to be used in diagnostic or clinical procedures.
Batch No.:	See product labels
Antigen:	The immunogen and standard used in this kit is recombinant human G-SCF amino acids 31-204 produced in and purified from an E.coli expression system.
Other Names:	Pluripoietin; G-CSF;
Accession:	P09919 CSF3_HUMAN;
Specificity:	Human G-CSF
Storage:	Store at 2-8C
Kit components:	The ELISA kit box contains 2 x 96-well pre-coated strip plates, protein standards, detection reagents, substrate buffer and detailed protocols.
Range:	31.2 pg/mL - 2,000 pg/mL
Sensitivity:	< 4 pg/ml
Kit protocol:	Please refer to our online product listing for current protocol/MSDS versions.

IX: Typical Standard Curve (for reference only, not to be used for actual data)

Concentration pg/ml	0.0	31.2	62.5	125	250	500	1000	2000
O.D.	0.012	0.183	0.329	0.558	0.860	1.252	1.642	1.922



This standard curve is for demonstration purposes only. A standard curve should be generated for each assay.

FOR RESEARCH USE ONLY