

Anti-Yellow Fluorescent Protein (Mouse Monoclonal)

Catalogue No.:	M-1320-100
Batch No.:	See product label
Unit size:	100 ug
Antigen:	Green Fluorescent Protein (GFP) from the jellyfish <i>Aequorea victoria</i> N-Terminal peptide-KLH conjugates.
Other Names:	YFP;
Produced in:	Mouse
Purity:	Affinity purified
Applications:	Western Blotting (WB), ELISA, Dot Blot, Immunoprecipitation, Immunostaining. Western Blotting suggested dilution at 1:1000. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
Specificity:	Recognizes native and denatured forms of Green Fluorescent Protein (GFP) and its variants: Enhanced Green Fluorescent Protein (EGFP), Yellow Fluorescent Protein (YFP), Enhanced Yellow Fluorescent Protein (EYFP) and Cyan Fluorescent Protein (CFP).
Form:	Liquid. PBS, pH 7.4 with 0.02% sodium azide
Storage:	Stable for 1 year at -20degC from the date of receipt unopened. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. After opening maintain in undiluted aliquots for up to 6 months. Avoid repeated freezing and thawing.
Expiry Date:	12 months after date of receipt, unopened.

FOR RESEARCH USE ONLY