

Chicken polyclonal antibody to acetyl Lysine: Affinity purified

Catalogue No.:	C-1516-100
Description:	Lysine acetylation of histones and non-histone proteins plays an important part in many cellular processes such as chromatin and nuclear signaling, transcription, gene silencing, cell cycle progression, apoptosis, differentiation, DNA replication and repair.
Batch No.:	See product label
Unit size:	100 uL
Antigen:	Acetylated Lysine OVA
Other Names:	Acetyl-lysine;
Produced in:	Chicken
Applications:	ELISA and WB. Suggested dilution of 1:1,000-1:5,000. Biosensis recommends that the optimal working dilution should be determined by the end user.
Specificity:	Binds to proteins with acetylated lysine residues.
Form:	Liquid. PBS with 0.02% Sodium Azide
Storage:	Short term storage at 2-8C for one week. At -20C as an undiluted liquid for up to 12 months.
Expiry Date:	12 months after purchase

FOR RESEARCH USE ONLY