

Chicken polyclonal antibody to Adiponectin (230-244): Affinity purified

Catalogue No.:	C-1524-100
Description:	Adiponectin is synthesized by adipocytes and is involved in the control of fat metabolism and insulin sensitivity, with direct anti-diabetic, anti-atherogenic and anti-inflammatory activities (ref: SWISSPROT).
Batch No.:	See product label
Unit size:	100 uL
Antigen:	A synthetic peptide from human Adiponectin (230-244 aa).
Sequence:	NDSTFTGFLLYHDTN
Other Names:	Adipocyte; C1q and collagen domain-containing protein; ADIPOQ; Adipocyte complement-related 30 kDa protein; ACRP30; ACDC; ACRP30; APM1; GBP28;
Accession:	Q15848 ADIPO_HUMAN;
Produced in:	Chicken
Applications:	WB. Suggested dilution of 1:2,000-1:5,000. Biosensis recommends that the optimal working dilution should be determined by the end user.
Specificity:	This antibody detects adiponectin protein at approx 30 kDa in human, mouse and rat.
Cross-reactivity:	Human, Mouse, Rat
Form:	Liquid. PBS with 0.02% Sodium Azide
Storage:	Short term storage at 2-8C for one week. At -20C as an undiluted liquid for up to 12 months.
Expiry Date:	12 months after purchase

FOR RESEARCH USE ONLY