

## Chicken polyclonal antibody to human Fatty acid-binding protein, adipocyte (97-111): Affinity purified

<b>Catalogue No.:</b>	C-1527-100
<b>Description:</b>	Fatty acid-binding protein, adipocyte (ALBP) is a lipid transport protein which binds long chain fatty acids and other hydrophobic ligands and delivers them to their receptors in the nucleus. ALBP is found in the cytoplasm and nucleus of adipocytes.
<b>Batch No.:</b>	See product label
<b>Unit size:</b>	100 uL
<b>Antigen:</b>	Human Fatty acid-binding protein, adipocyte peptide (97-111 aa).
<b>Sequence:</b>	KWDGKSTTIK RKRED
<b>Other Names:</b>	ALBP; Fatty Acid Binding protein; aP2; Adipocyte lipid-binding protein; Adipocyte-type fatty acid-binding protein; A-FABP; AFABP; Fatty acid-binding protein 4; FABP4;
<b>Accession:</b>	P15090 FABP4_HUMAN;
<b>Produced in:</b>	Chicken
<b>Applications:</b>	WB. Suggested dilution of 1:2,000-1:5,000. Biosensis recommends that the optimal working dilution should be determined by the end user.
<b>Specificity:</b>	This antibody detects human Fatty acid-binding protein, adipocyte at approx 16 kDa.
<b>Cross-reactivity:</b>	Human
<b>Form:</b>	Liquid. PBS with 0.02% Sodium Azide
<b>Storage:</b>	Short term storage at 2-8C for one week. At -20C as an undiluted liquid for up to 12 months.
<b>Expiry Date:</b>	12 months after purchase

---

FOR RESEARCH USE ONLY

---