

Human Plasminogen activator inhibitor 1 ELISA Kit (2 plates)

Catalogue No.: BEK-2142-2P

Description: The human Plasminogen activator inhibitor 1 (PAI-1) Kit is a sandwich ELISA. The capture antibody is a polyclonal human PAI-1 antibody pre-coated onto the 96-well strip plates provided in the kit. Human test samples and standards of known PAI-1 concentration are added to these wells and allowed to complex with the bound PAI-1 antibody. A biotinylated human PAI-1 polyclonal antibody is then added. This detection antibody binds to the antigen thus completing the sandwich. After washing, an enzyme Avidin-Biotin-Peroxidase complex (ABC) is added which binds to the second antibody. The peroxidase substrate TMB is added to induce a coloured reaction product. The intensity of this coloured product is directly proportional to the concentration of PAI-1 present in the samples. The purpose of this kit is the in-vitro quantitative determination of human PAI-1 in samples such as sera, plasma, tissue lysates and cell culture supernates. This kit has been configured for research use only and is not to be used in diagnostic or clinical procedures.

Batch No.: See product labels

Other Names: Endothelial plasminogen activator inhibitor; Serpin E1; SERPINE1; PAI1; PLANH1; PAI-1;

Accession: P05121 PAI1_HUMAN;

Specificity: Human Plasminogen activator inhibitor 1

Storage: Store at 2-8C

Kit components: The ELISA kit box contains 2 x 96-well pre-coated strip plates, protein standards, detection reagents, substrate buffer and detailed protocols.


Range: 156 pg/mL - 10,000 pg/mL

Sensitivity: < 10 pg/ml

Kit protocol: Please refer to our online product listing for current protocol/MSDS versions.

Concentration pg/ml	0.0	0.1	156	312	625	1250	2500	5000	10,000
O.D.		0.023	0.111	0.205	0.380	0.725	1.192	1.594	2.061

This standard curve is for demonstration purposes only. A standard curve should be generated for each assay.



FOR RESEARCH USE ONLY