

Mouse CD30 ligand ELISA Kit (2 plates)

Catalogue No.: BEK-2164-2P

Description: The mouse CD30 ligand (CD30L) Kit is a sandwich ELISA. The capture antibody is a monoclonal mouse CD30L antibody pre-coated onto the 96-well strip plates provided in the kit. Mouse test samples and standards of known CD30L concentration are added to these wells and allowed to complex with the bound CD30L antibody. A biotinylated mouse CD30L polyclonal antibody is then added. This detection antibody binds to the antigen thus completing the sandwich. After washing, an enzyme Avidin-Biotin-Peroxidase complex (ABC) is added which binds to the second antibody. The peroxidase substrate TMB is added to induce a coloured reaction product. The intensity of this coloured product is directly proportional to the concentration of CD30L present in the samples. The purpose of this kit is the in-vitro quantitative determination of mouse CD30L in samples such as sera, plasma, tissue lysates and cell culture supernates. This kit has been configured for research use only and is not to be used in diagnostic or clinical procedures.

Batch No.: See product labels

Other Names: Tumor necrosis factor ligand superfamily member 8; CD153; Tnfsf8;

Accession: P32972 TNFL8_MOUSE;

Specificity: Mouse CD30L

Storage: Store at 2-8C

Kit components: The ELISA kit box contains 1 x 96-well pre-coated strip plate, protein standards, detection reagents, substrate buffer and detailed protocols.

Range: 15.6 pg/mL - 1,000 pg/mL


Sensitivity: < 1 pg/ml

Kit protocol: Please refer to our online product listing for current protocol/MSDS versions.

Typical Mouse CD30L ELISA Kit Standard Curve

(TMB reaction incubate at 37°C for 16 min)

Concentration	0.1pg/ml	15.6pg/ml	31.3pg/ml	62.5pg/ml	125pg/ml	250pg/ml	500pg/ml	1000pg/ml
O.D.	0.064	0.102	0.143	0.232	0.393	0.699	1.301	2.369



This standard curve is for demonstration purposes only. A standard curve should be generated for each assay.

FOR RESEARCH USE ONLY