


Mouse monoclonal antibody to human AHNAK nucleoprotein [3G7]: IgG

Catalogue No.:	M-824-100
Batch No.:	See product label
Unit size:	100 ug
Antigen:	Partial recombinant human AHNAK protein (1-101) with a GST tag.
Clone:	3G7
Other Names:	Desmoyokin; AHNAK
Accession:	A1A586_HUMAN
Produced in:	Mouse
Purity:	IgG purified monoclonal antibody.
Applications:	This antibody is recommended for use in WB , direct ELISA and sandwich ELISA. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
Specificity:	The specificity of this antibody was demonstrated by WB against the recombinant protein.
Cross-reactivity:	This antibody has not yet been tested against other species.
Form:	Lyophilised from PBS pH 7.2, 150mM NaCl.
Reconstitution:	Reconstitute with 100 uL of sterile water. Centrifuge to remove any insoluble material.
Storage:	After reconstitution keep aliquots at -20C for a higher stability, and at 2-8C with an appropriate antibacterial agent. Glycerol (1:1) may be added for an additional stability. Avoid repetitive freeze/thaw cycles.
Expiry Date:	12 months upon purchase.



Detection of the partial AHNAK recombinant protein in transfected 293T cell lysates with the mouse monoclonal antibody to human AHNAK nucleoprotein [3G7], catalogue number M-800-100. This antibody detects a 16kDa band (lane 1) corresponding to the expressed recombinant partial human AHNAK protein. The antibody is not-reactive with non-transfected lysate (lane 2).

FOR RESEARCH USE ONLY