

Mouse monoclonal antibody to human protein AC1 [3F12]: IgG

Catalogue No.:	M-808-100
Description:	The unprocessed precursor of protein AC1 has a length of 93 amino acids and an estimated molecular weight of 10.5KDa. This protein is expressed in neuroblastoma but its function is undetermined.
Batch No.:	See product label
Unit size:	100 ug
Antigen:	Recombinant protein of human protein AC1 (amino acids 1 to 93) with a GST tag.
Clone:	3F12
Other Names:	Uncharacterised protein c4orf6
Accession:	CD006_HUMAN
Produced in:	Mouse
Purity:	Protein G purified immunoglobulin
Applications:	This antibody is recommended for ELISA. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
Specificity:	The specificity of this antibody was shown by ELISA using the recombinant human protein AC1 with GST tag.
Cross-reactivity:	Other species have not yet been tested.
Form:	Lyophilised from PBS pH 7.2, 150mM NaCl.
Reconstitution:	Reconstitute in 100 uL of sterile water. Centrifuge to remove any insoluble material.
Storage:	After reconstitution keep aliquots at -20C for higher stability or at 2-8C with an appropriate antibacterial agent. Glycerol (1:1) may be added for additional stability. Avoid repetitive freeze/thaw cycles.
Expiry Date:	12 months after purchase

FOR RESEARCH USE ONLY