

## Mouse monoclonal antibody to PSTAIR [IL-16]: IgG

<b>Catalogue No.:</b>	M-1219-100
<b>Description:</b>	Cyclin dependent kinases (CDK) contain an evolutionary conserved PSTAIR sequence (EGVPSTAIRESLLKE) which distinguishes them from other protein kinases.
<b>Batch No.:</b>	See product label
<b>Unit size:</b>	100 ug
<b>Antigen:</b>	Synthetic 16 amino acid oligopeptide containing the PSTAIR sequence (EGVPSTAIRESLLKE) conjugated to BSA
<b>Other Names:</b>	EGVPSTAIRESLLKE antibody; PSTAIRE;
<b>Produced in:</b>	Mouse
<b>Purity:</b>	IgG
<b>Applications:</b>	Western Blotting (WB). A concentration of 0.25-0.50 ug/mL is recommended.
<b>Specificity:</b>	Human; mouse;
<b>Cross-reactivity:</b>	The specificity of this antibody has been confirmed by WB against the antigen.
<b>Form:</b>	Lyophilized from 1.2% sodium acetate, 2mg BSA, 0.01mg NaN <sub>3</sub>
<b>Reconstitution:</b>	Reconstitute in 1 mL of PBS (pH 7.4) to achieve an antibody concentration of 100 ug/mL. Centrifuge to remove any insoluble material.
<b>Storage:</b>	At least 12 months after purchase at 2-8C (lyophilized formulations). After reconstitution, aliquot and store at -20C for a higher stability. Avoid freeze-thaw cycles.
<b>Expiry Date:</b>	12 months after purchase.

---

FOR RESEARCH USE ONLY