


## Mouse polyclonal antibody to human NPIP: whole serum

<b>Catalogue No.:</b>	Mo-939-50
<b>Description:</b>	The deduced 350-amino acid NPIP protein is membrane bound. Immunofluorescence microscopy demonstrated localization of NPIP to the nuclear membrane and colocalization with nuclear pore complex protein nucleoporin p62. NPIP is widely expressed and is thought to associate with the nuclear pore complex.
<b>Batch No.:</b>	See product label
<b>Unit size:</b>	50 $\mu$ L
<b>Antigen:</b>	Partial recombinant human NPIP (97-196) with a GST tag.
<b>Other Names:</b>	Nuclear pore complex-interacting protein; NPIP
<b>Accession:</b>	NPIP_HUMAN
<b>Produced in:</b>	Mouse
<b>Purity:</b>	Whole serum
<b>Applications:</b>	This antibody is recommended for WB and ELISA. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
<b>Specificity:</b>	Specificity has been confirmed by WB and direct ELISA against the antigen.
<b>Cross-reactivity:</b>	Human. Other species have not been tested.
<b>Form:</b>	Liquid. 50% glycerol
<b>Storage:</b>	Keep aliquots at -20C for higher stability or at 2-8C with an appropriate antibacterial agent. Avoid repetitive freeze/thaw cycles.
<b>Expiry Date:</b>	12 months after purchase



Western blot detection of GST tagged recombinant human NPIP. Note the GST tag alone is 26 kDa.

FOR RESEARCH USE ONLY