

## Rabbit antibody to Noxa (1-16): whole serum

<b>Catalogue No.:</b>	R-123-100
<b>Description:</b>	The Bcl-2 family of proteins which regulate apoptosis share identical sequences called Bcl-2 Homology domains (BH1-4). The BH3 proteins, including BID, NOXA, PUMA, BIK, BIM and BAD are all pro-apoptotic and share sequence identity within the amphipathic alpha-helical BH3 region, which is essential for their apoptotic function. NOXA is highly expressed in adult T-cell leukemia cell line.
<b>Batch No.:</b>	See product label
<b>Unit size:</b>	100 uL
<b>Antigen:</b>	A synthetic peptide (MPGRKARRNA PVNPTR) as part of mouse Noxa (aa: 1-16) conjugated to diphtheria toxoid
<b>Other Names:</b>	PMAIP1; phorbol-12-myristate-13-acetate-induced protein 1; adult T cell leukemia-derived PMA-responsive; Immediate-early-response protein APR; PMA-induced protein 1; Pmaip1; Noxa
<b>Accession:</b>	NOXA_MOUSE
<b>Produced in:</b>	Rabbit
<b>Purity:</b>	Whole serum
<b>Applications:</b>	WB. A dilution of 1:1000 to 1:2000 is recommended for this application. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
<b>Specificity:</b>	Western blot analysis of cells infected with Noxa adenoviruses and BAF indicates a high level of specificity for this antiserum.
<b>Cross-reactivity:</b>	This antiserum cross-reacts with mouse. Not yet tested in other species.
<b>Form:</b>	Lyophilised
<b>Reconstitution:</b>	Reconstitute in 100 uL of sterile water. Centrifuge to remove any insoluble material.
<b>Storage:</b>	After reconstitution keep aliquots at -20C for a higher stability, and at 2-8C with an appropriate antibacterial agent. Glycerol (1:1) may be added for an additional stability. Avoid repetitive freeze/thaw cycles.
<b>References:</b>	<ol style="list-style-type: none"><li>1. Jansson, A.K., et al., Oncogene 22(30):4675-4678 (2003).</li><li>2. Hijikata, M., et al., J. Virol. 64(10):4632-4639 (1990).</li></ol>

---

FOR RESEARCH USE ONLY