

Rabbit polyclonal antibody to human Matrix metalloproteinase-8 (103-115): Affinity purified

Catalogue No.:	R-1125-100
Description:	The matrix metalloproteinases (MMPs) are a large family of zinc endopeptidases. All MMPs are synthesized as inactive proenzymes. The activation of these proenzymes is a critical step that leads to degradation of extracellular matrix components such as fibronectin and collagen type III. MMP8 is known to degrade fibrillar type I, II, and III collagens.
Batch No.:	See product label
Unit size:	100 ug
Antigen:	A synthetic peptide (PGNPKWERTNLTY) corresponding to a region (103-115) from human Matrix metalloproteinase-8. To enhance the immunological response, this peptide was coupled to carrier protein BSA.
Other Names:	Neutrophil collagenase; EC 3.4.24.34; Matrix metalloproteinase-8; MMP-8; PMNL collagenase; PMNL-CL; MMP8; CLG1;
Accession:	P22894 MMP8_HUMAN;
Produced in:	Rabbit
Purity:	Affinity purified on antigen column
Applications:	Immunohistochemistry (IHC) and Western Blotting (WB). A concentration of 0.75 ug/mL is recommended for WB. Human MMP8 (precursor) has a predicted length of 467 residues and MW of 53 kDa. A concentration of 1.0-2.0 ug/mL is recommended to detect the protein in formalin fixed and paraffin embedded tissues. Biosensis recommends optimal dilutions/concentrations should be determined by the end user.
Specificity:	The specificity of this antibody has been confirmed by WB and IHC against the antigen.
Cross-reactivity:	Human; rat; predicted to react with mouse due to sequence homology;
Form:	Lyophilised with 5mg BSA, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Thimerosal, 0.05mg NaN ₃
Reconstitution:	Reconstitute in 100 uL of sterile distilled water to achieve an antibody concentration of 1 mg/mL. Centrifuge to remove any insoluble material.
Storage:	At least 12 months after purchase at 2-8C (lyophilized formulations). After reconstitution, aliquot and store at -20C for a higher stability. Avoid freeze-thaw cycles
Expiry Date:	12 months after purchase

FOR RESEARCH USE ONLY