


## Rat Interleukin 2 ELISA Kit (2 plates)

- Catalogue No.:** BEK-2168-2P
- Description:** The rat Interleukin 2 (IL-2) Kit is a sandwich ELISA. The capture antibody is a monoclonal rat IL-2 antibody pre-coated onto the 96-well strip plates provided in the kit. Rat test samples and standards of known IL-2 concentration are added to these wells and allowed to complex with the bound IL-2 antibody. A biotinylated rat IL-2 monoclonal antibody is then added. This detection antibody binds to the antigen thus completing the sandwich. After washing, an enzyme Avidin-Biotin-Peroxidase complex (ABC) is added which binds to the second antibody. The peroxidase substrate TMB is added to induce a coloured reaction product. The intensity of this coloured product is directly proportional to the concentration of IL-2 present in the samples. The purpose of this kit is the in-vitro quantitative determination of rat IL-2 in samples such as sera, plasma, tissue lysates and cell culture supernates. This kit has been configured for research use only and is not to be used in diagnostic or clinical procedures.
- Batch No.:** See product labels
- Other Names:** Interleukin-2; IL-2; T-cell growth factor; TCGF; IL2; IL-2;
- Accession:** P17108 IL2\_RAT;
- Storage:** Store at 2-8C
- Kit components:** The ELISA kit box contains 1 x 96-well pre-coated strip plate, protein standards, detection reagents, substrate buffer and detailed protocols.
- Range:** 15.6 pg/mL - 1,000 pg/mL
- Sensitivity:** < 1 pg/ml
- Kit protocol:** Please refer to our online product listing for current protocol/MSDS versions.

Concentration pg/ml	0.0 pg/ml	31.2 pg/ml	62.5 pg/ml	125 pg/ml	250 pg/ml	500 pg/ml	1000 pg/ml	2000 pg/ml
O.D.	0.133	0.190	0.226	0.342	0.528	0.941	1.415	2.177

This standard curve is for demonstration purposes only. A standard curve should be generated for each assay.



FOR RESEARCH USE ONLY